


**FASTER
EASIER
BETTER!!!**



**Presenting Easy To
Use E-Z Prep QBC
Tubes!**



EZ Prep Pipetter

Professional Medical has discovered an easier way to perform cell counts on your Autoread and QBC-V. With the EZ Prep the float is already in the tube and the stopper is already partially seated. You simply fill the tube as you normally would, using your EZ Prep Pipetter.

After the tube is filled, simply seat the stopper, roll the tube gently and insert it into the QBC centrifuge. Your hands-on time is significantly reduced, allowing you to increase productivity. And since the EZ Prep packaging is about 1/3 its former size, you'll save on storage space.

With Professional Medical you can purchase as few as 10 QBC EZ Prep Tubes at a time!!

10 QBC EZ Prep Tubes Only \$42.50

Professional Medical Sales

520 N Parramore Avenue • Orlando, FL 32801

407-843-9797 • Toll Free 888-VET-STUF that's

888-838-7883

FLOAT

*Already in
the Tube.*

STOPPER

*Partially
Seated*

